# Information for Patients

## PERLE Mammary Implants

<b>EN</b>	English
<b>FR</b>	French
<b>DE</b>	German
<b>IT</b>	Italian
<b>ES</b>	Spanish
<b>PT</b>	Portuguese
<b>DA</b>	Danish
<b>NL</b>	Dutch
<b>EL</b>	Greek
<b>LV</b>	Latvian
<b>SK</b>	Slovak
<b>RO</b>	Romanian
<b>TR</b>	Turkish
<b>RU</b>	Russian
<b>CS</b>	Czech
<b>ET</b>	Estonian
<b>HU</b>	Hungarian
<b>PL</b>	Polish
<b>SL</b>	Slovenian
<b>SV</b>	Swedish
<b>LT</b>	Lithuanian
<b>BG</b>	Bulgarian
<b>HR</b>	Croatian
<b>FI</b>	Finnish
<b>NO</b>	Norwegian
<b>IS</b>	Icelandic

## EN - Information for patients

### Implant Card Purpose

The implant card filled by the health care professional should be provided to you after the surgery. The implant card must be kept on your person at all times and enables you:

- To identify the implanted devices and to get access to other information related to the implanted device (e.g. via EUDAMED, and other websites).
- To identify yourselves as persons requiring special care in relevant situations e.g. medical examinations, security checks.
- To inform emergency clinical staff or first responder about special care/needs for you in case of emergency situations.

The International Implant Card is linked with the GCA Comfort™ Plus Warranty.

This card applies to all GC Aesthetics® Silicone Gel Filled Breast Implants (certified under MDR 2017/745) implanted by an appropriately qualified and licensed Physician in accordance with accepted surgical procedures and techniques, and in strict accordance with the GCA® "Instruction for Use" supplied with the product.

### Patient leaflet Purpose

The purpose of this leaflet is to provide additional information you should be aware in terms of:

1. Information about the **expected lifetime of the device** and any necessary follow up.
2. Any information to **ensure safe use of the device by the patient**, including information on the materials and substances to which you can be exposed.
3. **Warnings, precautions or measures** to be taken by the patient or a healthcare professional with regard to interference with any possible external influences, medical examinations or environmental conditions.

The updated version of this leaflet is available on our information website for patients, by the link detailed on the implant card:

<https://www.gcaesthetics.com/patient>

#### 1. Product Lifetime

Breast implants are not lifetime devices. Be aware that they may need to be removed or replaced at some future date which can lead to revision surgery. Many surgeons suggest replacement after 10-20 years. But unless there is a problem, this may not be necessary. But unless there is a problem, this may not be necessary.

#### 2. Necessary Follow-Up / Information to Ensure Safety

##### Postoperative Care

Your surgeon will provide advice on how to take care of yourself after you have breast implant surgery. These may include you wearing a supportive bra 24 hours a day or changing the adhesive bandages when necessary. Your surgeon may also prescribe prophylactic antibiotics. These and other measures are at your surgeon's opinion and are based on their personal experience and overall treatment plan.

Not following the surgeon's recommendations may lead to complications.

You should consult your surgeon if you suspect any complications. In the case of trauma to the area, please speak to your surgeon. Trauma can include extreme massaging of the breast regions, by some sporting activity or by using seat belts. Contact sports are not recommended. Tanning at salons or sunbathing can influence scar healing. Direct scar exposure to the sun should be avoided for about six months after surgery. It is recommended that you consult a physician or pharmacist before the use of topical medicines (e.g., steroids) in the breast area.

If any surgery of the breast area is scheduled, you should inform the physician or the surgeon of the presence of an implant and carry your implant card.

The postoperative care recommendations presented within this section of the patient brochure does not include all possible options. Your surgeon should provide postoperative care details.

##### After Surgery

Once breast augmentation surgery has taken place, a follow-up appointment with a surgeon will be scheduled. This is routine and will be to check on your general health after the operation.

Following this, contact with a surgeon should be limited if you continue to feel fit and healthy.

However, you should always remember to consult a surgeon, physician or pharmacist before the use of topical medicines (e.g. steroids) in the breast area. You should also inform them of the presence of implants if you require any surgery in the future.

##### Self-Examination and Auto palpation

It's recommended that adult women of all ages perform a breast self-exam at least once a month. A regular self-exam will help you to be familiar with how your breasts look and feel. You can contact a healthcare professional if you notice any changes.

One of the ways you can perform a breast self-exam is in front of a mirror by looking for any changes in contour, any swelling, or dimpling of the skin, or changes in the nipples. You should do this with your arms at your sides and your arms high overhead. Next, rest your palms on your hips and press firmly to flex your chest muscles. Look for any dimpling, puckering, or changes, particularly on one side.

##### Breast and Nipple Piercing Procedures

Women with breast implants who undergo body piercing procedures to the breast region must consider the possibility that an infection could develop following this procedure. Treatment including antibiotics, possible removal of the implant, or additional surgery may be necessary.

##### Materials

These medical devices are manufactured from medical grade silicone materials which are suitable for long term implantation. Master Files have been filed with the U.S FDA.

The following three tables provide quantitative and qualitative information on the materials and substances which patients can be exposed based on chemical characterization of representative devices. The materials and device have been subjected to biocompatibility testing and evaluation and risk assessments to demonstrate their biological safety. However, individual responses to chemicals may vary, and all reactions cannot be predicted.

**Table 1. PERLE Mammary Implant device materials**

Device Materials	Implant component
Silicone dispersions	Shell
Silicone dispersions	Closure patch
Silicone adhesive	shell
Silicone gel	Gel fill

**Table 2. chemicals released from the PERLE Mammary Implant**

Compounds	Whole device (µg/g or ppm)	Compounds	Whole device (µg/g or ppm)
<b>Volatiles<sup>1</sup> – chemicals that can be released by the mammary implant as a gas</b>			
Benzene	Not Detected	Ethylbenzene	Not Detected
Bromobenzene	Not Detected	Hexachlorobutadiene	Not Detected
Bromodichloromethane	Not Detected	Isopropylbenzene	Not Detected
Bromoform	Not Detected	p-Isopropyltoluene	Not Detected
Butylbenzene (N-/tert-/sec)	Not Detected	N-propylbenzene	Not Detected
Carbon tetrachloride	Not Detected	Naphthalene	Not Detected
Chlorobenzene	Not Detected	Styrene	Not Detected
Chloroform	Not Detected	1,1,1,2-Tetrachloroethane	Not Detected
2,4-Chlorotoluene	Not Detected	Tetrachloroethene	Not Detected
Dibromochloromethane	Not Detected	Trichloroethene	Not Detected
1,2-Dibromo-3-chloropropane	Not Detected	Trichloroethane (1,1,1-/1,1,2-)	Not Detected
Dichlorobenzene(1,2-/1,3-/trans-1,2-)	Not Detected	1,2,3-Trichloropropane	Not Detected
Dichloroethane (1,1-/1,3-)	Not Detected	Trimethylbenzene(1,2,4-/1,3,5-)	Not Detected
Dichloroethene (1,1-/ cis-1,2-/ trans-1,2-)	Not Detected	Toulene	Not Detected
Dichloropropane (1,2-/1,3-)	Not Detected	Trichlorobenzenes (1,2,4-/ 1,2,3-)	Not Detected
1,1-Dichloropropene	Not Detected	o-Xylene	Not Detected
Ethanol, 2-(trimethylsilyl)-	0.036 max.	m- & p-Xylene	Not Detected
Total Volatiles	0.036		
<b>Non-volatiles</b>			
<b>chemicals that can be released by breast implants following soaking in water</b>			
Chlorinated aromatic compound	0.08	mono chlorinated aromatic compound	0.08
<b>chemicals that can be released by breast implants following soaking in solvent</b>			
mono chlorinated aromatic compound	0.16	C-9 dialkyl adipate	0.05
Irgafos 168 oxide fragment: one ether bond hydrolyzed	0.08	irgafos 168 oxide	0.05
polyoxygenated unsaturated compound	0.08	C-8 dialkyl benzene dicarboxylate <sup>2</sup>	2.51
C-8 dialkyl benzene dicarboxylate	1.65		
Total non-volatiles	4.74 Max		
<b>Extractables</b>			
<b>chemicals that can be released by breast implants following soaking in water<sup>3</sup></b>			
Cyclic siloxanes (D4, D5, D6..)	Not Detected	Linear siloxanes (L3, L4, L5...)	Not Detected
<b>chemicals that can be released by breast implants following soaking in solvent</b>			
Cyclic siloxanes (D4, D5, D6..)	1009.02 max	Linear siloxanes (L3, L4, L5...)	8.1 Max
Total extractables	≤1017.12		
<b>Other compounds of interest</b>			
2,4-Dichlorobenzoyl peroxide	0.1 Max		
<sup>1</sup> Not Detected means that the level of the individual volatile was below the quantification limit of the test method. The quantification limit for all compounds was 0.006 µg/g except for m/p-xylene and 2/4-chlorotoluene, which have quantification limit of 0.012µg/g			
<sup>2</sup> detected using dichloromethane (DCM) solvent extraction			
<sup>3</sup> Not Detected means that the level of individual extractable was below 6.12µg/g, the quantitation limit of the test method.			

**Table 3. Heavy Metals found in the PERLE Mammary Implants**

Heavy metals	Concentration (ppm)	Heavy Metals	Concentration (ppm)
Arsenic	Not Detected	Chromium	Not Detected
Cadmium	Not Detected	Copper	Not Detected
Lead	Not Detected	Tin	Not Detected
Mercury	10 max.	Nickel	Not Detected
Vanadium	Not Detected	Platinum	Not Detected
Molybdenum	Not Detected	Zinc	0.002 max.
Selenium	Not Detected	Magnesium	Not Detected
Cobalt	Not Detected	Manganese	Not Detected
Antimony	Not Detected	Thallium	Not Detected
Tungsten	Not Detected		
Not Detected means that the level of the individual element was below 0.0012ppm, the quantification limit of the test method.			

Based on certification from suppliers the devices are manufactured from materials which:

- Do not contain Carcinogens, Mutagens and Reproductive Toxic Substances (CMR) or Endocrine- Disruptive (ED) substances above regulated threshold levels requiring specific labelling
- Do not intentionally contain Phthalates
- Do not incorporate materials of biological origin

### 3. Warnings, Precautions to be Taken

#### Mammography

Women who have breast implants should continue to talk to their physician to carry out normal checks to detect breast cancer. The implant can hinder the detection of early breast cancer through mammography by concealing some underlying breast tissue. The implant can also compress overlying tissue which can 'hide' suspicious lesions in the breast.

Please make sure you share all information you know about your breast implants with your healthcare provider so that they can adapt the mammography pressure accordingly. You should request radiologists who are experienced with the most current techniques and equipment for imaging breasts with implants.

You should inform them of the presence, type and placement of implants and the orientation marks on the anatomical implant. The presence of a breast implant may make screening mammography more difficult and may lead to additional mammography views.

#### Detection of Breast Cancer

The presence of breast implants could delay the detection of breast cancer by self-exam. For this reason, patients must be informed that they need to consult a surgeon or doctor for appropriate medical monitoring. Patients should also get regular breast cancer screening.

Possible Interferences of the Implant in MR Environment / Security Scanners.

Please note that Eurosilicone's Gel-filled Breast Implants are all made from medical implant grade silicone materials which can be used with MRI Scan. No precaution to be taken.

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© 2022 GC Aesthetics.® 904767.
GC Aesthetics® is a registered trademark
Global Consolidated Aesthetics Limited, registered in Ireland with number 450181.
of Global Consolidated Aesthetics Limited.